

## PCR (Saliva) Examination ※Collection Procedure

Please refrain from eating, drinking, smoking, brushing your teeth, and rinsing your mouth 30 minutes before collection.

- A. Remove the cap of the provided container (1), then , attach the funnel (2).
- B. Press the tip of your mouth on the top of the funnel, and slowly pour saliva into the container (3).
- C. Gently pour the saliva into the funnel. (Please spit out at least 2.0ml of saliva.)  
When the prescribed amount of saliva has been collected, remove the funnel and reattach the cap.
- D. Put this container in the submission bag, and close the double zipper bags tightly (4).
- E. Then your specimen bag is dropped into the“PCR BOX” located back entrance of clinic (5).  
(cf,MAP , back of this page).
- F. After examination,this funnel and plastic bag, please bring with you, when you return (6).

(1) Container cap



(2) Funnel




(3) Collection of Saliva



(4) Double zipper bags



(5)



(6)



